

12611 Non-Vented

Full Face Mask

INTENDED OF USE :

This Mask is intended to be used under the specific direction of a physician. It is used in conjunction with Non Invasive Ventilators in hospital and Institutional environments.

STRUCTURE AND MATERIALS :

Headgear : CR LATEX, SBR and Nylon Yarn
Buckle : PC
Connector : PC
O ring : Silicone
Frame : PC
Cushion : Silicone

LATEX FREE

CONNECTIONS COMPATIBILITY : ISO 5356-1:2015

SINGLE PATIENT USE

STERILITY : Non-sterilized

AUTOCLAVABLE : Autoclavable, up to 121°C/250°F

PRODUCT VALIDITY : 5 years after the manufacture date

CONSERVATION PROCEDURE :

Suggested stock temperature : 15°C ~ 25°C
Storage R.H.: <85%
To be stocked in a dry and clean place
Avoid Direct Sunray

