


DualSpa

Nasal Mask



- A** New 3D mask cushion provides an excellent seal without tight headgear; the face pressure sores and discomfort won't happen.
- B** Just a simple action to attach and detach the cushion from frame.
- C** **FinTek** >> Innovative adjustable forehead support provides stability with easy one - handed adjustment and allows the mask fit most of patients.
- D** Flexible, one-piece forehead pad evenly conforms to the forehead contour, so the pad stays in place and you won't feel a thing.
- E** 360° rotational ball in elbow.