


820

Full Face Mask Nasal Mask



- A** Medical grade liquid silicone full face cushion for greater comfort special 3D design for better sealing user's face but don't need to tighten the headgear too hardly.
- B** Quick assembling design offer user more convenience when taking off / assembling the cushion.
- C** Silence type vent holes provides lower noise and the maximum CO₂ washout, vented air is spread as cloud so that does not disturb patient or companion.
- D** Sampling port allows measuring pressure or ETCO₂.
- E** Anti-asphyxia valve, automatically open when the generator pressure fails, enabling patients to breathe by themselves. (Full Face Mask)
- F** Flexible, One-piece forehead pad evenly conforms to the forehead contour, so the pad stays in place and you won't feel a thing.
- G** Easy to unhook or reattach the upper headgear by threading it through the loop.